# **Rigel Uni-Therm**

Electrosurgical Analyser



Technical Data Sheet

Innovating Together





# Innovating Together

# Introducing the new high performance electrosurgical analyser with the emphasis on speed and ease of testing all modern day electrosurgical devices.

The new high power Rigel Uni-Therm accurately measures the performance of electrosurgical generators. Measurements include; high frequency leakage, high power current and power distribution and patient return plate alarm testing.

Boasting a full colour screen and providing easy navigation, the Rigel Uni-Therm is easy to operate straight out of the box. The large colour screen provides step by step guidance to ensure the correct connection to the device under test.


A high power, high resolution load bank enables performance testing up to 6A RMS with a duty cycle of 100% up to 60 seconds. This makes the Uni-Therm an extremely versatile test instrument for calibrating and performance testing of all high power electrosurgical instruments on the market. The new Rigel Uni-Therm incorporates a range of custom settings that include a variety of test sequences and power distribution curves which make testing fast, easy and accurate.

Utilising the large full colour display, built-in memory and data management facilities, the Rigel Uni-Therm is the first fully stand-alone diathermy tester that can offer true user friendly operation and operator safety. The large array of internal resistors, ranging from 0-5100 $\Omega$  in  $5\Omega$  steps provide not only the most accurate and detailed power curves, the Rigel Uni-Therm also advocates safe working practise by providing all necessary resistors within the enclosure.

The return electrode monitoring test offers a unique isolated and electronically driven potentiometer to provide accurate and fast testing of alarms.

In addition, the Bluetooth-enabled technology allows the wireless connectivity of the tester with PCs and other accessory equipment for fast and convenient downloading of performance data and the uploading of manufacturers' power curves and manufacturer specific test programs. Other special features of the Rigel Uni-Therm include in-built graphics to display connection diagrams, program settings, power curves, alarm settings, data management and data analyses.

## **Custom Test Settings**




The unique colour user-interface provides both graphic as well as numeric representation of the test data


Showing manufacturer specific curves and

## Key features

- Fully compliant with IEC 60601-2-2 One instrument for full compliance testing offering peace of mind
- Accurate and safe Utilising full 10kV isolation on all measuring systems
- High power load bank Measure up to 6 A RMS with dutycycle up to 100% for 60 seconds
- High frequency leakage
   Easy to connect with onscreen help for each configuration
- Power distribution curves
  Variable load with full 10kV isoltion from 0 to 5100Ω in 5Ω steps Accurate, fast, and flexible
- Remote electrode monitoring testing Using electronic potentiometer range upto 500Ω in 1Ω steps with high and low alarms
- Stand-alone Not relying on PC or laptop, direct print facility via Bluetooth
- Automatic and manual test sequences For fast and effective (repeat) testing
- Stylish and rugged enclosure
  Small footprint ideal for in-situ testing
- Graphic colour user interface For fast and easy navigation and connection to DUT
- Future upgrade ready Download future upgrades from the web into your tester
- Prepared for PPM protocols Configured for automatic performance testing of a variety of parameters








allowable tolerances, the large colour screen provides instant performance indication.

Use the graphic representation to quickly determine whether the test data meets the set criteria.

	RUN P	OWER '	TEST	(DATA	VIEW)	
#	Ohms	W rms	mA rms	V rms	V peak	CF
1	400	265	813	327	498	1.5
2	500	268	733	367	551	1.5
3	600	267	668	403	594	1.4
4	700	269	620	437	650	1.4
5	800	268	579	464	683	1.4
6	900	267	546	491	728	1.4
7	1000	264	513	>500	765	1.4
GOTO GOTO						
GRAPH						


Use the numeric representation of the data to analyse parameters such as peak voltage, crest factor, power and current at each resistance point

## **Remote Electrode Monitoring**




Use the rotary encoder to set high and low alarm points as well as controlling the potentiometer, scaling up or down in manual or automatic mode. Then, simply capture the alarm using the dedicated fast key on the screen. Data will be stored onboard for future traceability.

## **HF Leakage**



Individual settings can be selected using the simple configuration template. The settings are automatically translated into a detailed scheme as shown below:



Detailed instructions on the leakage path and connections between each setting are provided in a clear and colour coded overview.

Colour coded connections correspond to the matching connection on the tester.

## **High Power Load Bank**



The newly designed load bank is capable of measuring currents up to 6 A RMS with a maximum duty cycle of up to 100% for 60 seconds. This ensures the Rigel Uni-Therm is compatible with all high power electrosurgical devices.

## **Applications**

- Calibration of ESU's
- To verify the correct performance of Electrosurgical devices during:
  - a. Routine testing
  - **b.** Type testing
  - **c.** Production line testing
- Unique training tool for testing on electrosurgical devices

## What comes in the box?

- Mains lead
- Instruction manual
- Application disc
- Bluetooth USB adaptor
- Calibration certificate



## Innovating Together

## **SPECIFICATIONS**

Power measurement True RMS value of applied waveform

Power rating 0 - 500W (RMS)

Duty cycle 100% up to 60 seconds

Load bank  $0 - 5115\Omega$  Resolution  $5\Omega$ 

Accuracy (1W + 5% of value)

Voltage (peak) 0 - 10kV (Peak) - Closed load only

Accuracy  $\pm 10\%$  of value.

Measurement is taken between the active and dispersive electrodes with closed load only.

 Voltage
 0 - 700V (RMS)

 Accuracy
 (2V + 2% of value)

 Current
 0 - 6000mA (RMS)

 Accuracy
 ±10mA or 2% of value

 Crest factor
 1.4 - 20 (Vpeak / V RMS)

The higher of the two peak voltage measurements is used for calculation

**RMS Bandwidth** 

Instrumentation only

30 Hz to 10 MHz (-3 dB)

With loads:

30 Hz to 2.5 MHz (-3 dB)

Variable loads  $0 - 5115\Omega$ , steps @  $5\Omega$  (1023 steps)

Accuracy 1% of F

Load array

Ceramic resistors (Non inductive)

Measurement delay

Foot switch delay selectable between

200 - 5000ms (10mSec resolution)

RF Leakage (High Frequency Leakage)

Active From active part to earth

Passive From plate - receptacle - to earth Load Variable see power measurement

Fixed 2 x  $200\Omega$ 

**Return Electrode Monitoring** 

Range 1-  $475\Omega$ , steps @  $1\Omega$  steps

Motor driven potentiometer

Alarm register High and low, manual confirmation

Ranging Manual or automatic

**Output Connectors** 

Remote foot switch control (CUT) 2 x 4mm - yellow, single relay contact

Remote foot switch control (COAG) 2 x 4mm - blue, single relay contact

High Frequency leakage Through 4mm sockets and power measurements

USB PC download / Future PC software
Oscilloscope output 5V/A, 100mA RF current minimum input



# Rigel's Med-eKit can include the following:



## 288 Electrical Safety Analyser

- Light, hand-held, battery operation
- Conform IEC 62353 / 60601/
   VDE 0751 / NFPA-99 / AS-NZS 3551
- Memory for up to 10,000 devices
- Bluetooth communication
- Full, semi automatic & manual tests



#### **BP-SIM NIBP Simulator**

- Light, hand-held, battery operation
- Adult & paediatric NIBP simulations
- Manufacturer specific O-curves
- Overpressure and leak test
- Memory for up to 10,000 devices



## SP-SIM SpO2 Simulator

- Light, hand-held, battery operation
- Tests probe and monitor all at once
- User configurable simulations
- Manufacturer R-curves
- Memory for up to 10,000 devices

## As well as:

- Patient Simulator
- Flow Analyser
- Defib Analyser
- PrinterBarcode Scanner
- Asset Management Software
- Non-Rigel Test Equipment



# Innovating Together

Isolation 10kV Isolation between measurement device

and enclosure

**Low Frequency Filter** 100 Hz filter to avoid low-frequency disturbance

or interference

Storage and Recall

Memory Approx 5,000 records
Output CSV and SSS format

**General Specifications** 

Dimensions 370 X 300 X 204mm

Weight 10 kg

Operating temperature  $15 \,^{\circ}\text{C}$  to  $35 \,^{\circ}\text{C}$  Storage temperature  $0 \,^{\circ}\text{C}$  to  $50 \,^{\circ}\text{C}$ 

Mains power 120/230 VAC +10%; 48 to 66 Hz, 35 VA

Part number: 398A910

## Also available

#### Other extras

- Med-eBase asset management software
- Barcode scanner
- Test lead set
- Protective travel case (by peli case)
- USB Keyboard

## From Rigel Medical

- Rigel 266 Plus Manual Safety Analyser
- Rigel 277 Plus Automatic Safety Analyser
- Rigel 288 Hand-held Safety Analyser
- Rigel 601 Checkbox
- Rigel UNI-SIM Vital Signs Simulator
- Rigel BP-SIM NIBP Simulator
- Rigel SP-SIM SpO2 Simulator
- Rigel 333 Patient Simulator
- Rigel Uni-Pulse Defibrillator Tester
- Rigel Multi-Flo Infusion Pump Tester
- Med-eBase Software Application

## From the Seaward Group

- Portable Appliance Testers
- IEC Lead Tester
- Insulation Resistance Testers
- RCD Testers
- Earth Loop Impedance Testers
- Installation Testers
- Multimeters
- Current Clamps
- Hipot Testers
- Earthbond Testers
- Micro Ohmmeters