



Innovating Together

# An introduction to Infusion Pump Testing

# Double your test capacity

Save time, save money



## Introducing the NEW Rigel Multi-Flo Infusion Pump Analyzer

The only multi-channel infusion pump tester that provides instant flow rate, volume and pressure measurements for unbeatable speed and accuracy.

FREE Guide  
to Infusion  
Pump Testing

For your free copy visit

[www.seaward-groupusa.com/multi-flo](http://www.seaward-groupusa.com/multi-flo)

Visit [www.seaward-groupusa.com/multi-flo](http://www.seaward-groupusa.com/multi-flo) to find out more  
Or call 813-886-2775

 **RIGEL**  
MEDICAL

Innovating Together

## Contents

<b>Foreword</b>	2	<b>7 Multi-Flo</b>	21
<b>1 Introduction</b>	2	7.1 Introduction	21
<b>2 Introduction and History of Infusion</b>	3	7.2 Testing Infusion Devices with the Multi-Flo	22
2.1 Infusion Circuit	5	7.2.1 Syringe Driver Connection	22
<b>3 Types of Infusion</b>	5	7.2.2 Volumetric Pump Connection	22
3.1 Elastomeric Pumps	6	7.2.3 Volume/Flow	22
3.2 Gravity Controllers	7	7.2.4 Occlusion	23
3.3 Syringe Drivers	7	7.2.5 PCA	24
3.4 Volumetric Pump	8	<b>8 Conclusion</b>	24
3.5 Patient Controlled Analgesia (PCA) Pump	8	<b>9 References</b>	25
3.6 Ambulatory Pumps	9		
<b>4 Infusion Characteristics</b>	9		
<b>5 Problems Associated with Infusion</b>	10		
5.1 Air Embolism	11		
5.2 Free Flow or Siphonage	11		
5.3 Occlusion	11		
5.4 Bolus	12		
5.5 Vein Closing After Infusion	12		
5.6 Infiltration or Tissueing	12		
<b>6 Testing Infusion Devices</b>	13		
6.1 Introduction	13		
6.2 Volume/Flow Rate	14		
6.2.1 Cylinders and Burette	14		
6.2.2 Weighing Scales	15		
6.2.3 Vernier Calipers/Dial Gauges	16		
6.2.4 Electronic Devices	16		
6.3 Occlusion and Alarm Pressure	17		
6.3.1 Analogue Pressure/Force Gauge	17		
6.3.2 Digital Pressure Meter	18		
6.3.3 Electronic/Automatic Devices	19		
6.4 PCA	19		
6.4.1 Electronic/Automatic Analyzers	19		
6.5 Trumpet Curve	19		

## Foreword

This booklet is written as a guideline for people involved in testing medical infusion devices. All reasonable care has been taken to ensure that the information, reference figures and data are accurate and have been taken from the latest versions of various standards, guidance notes and recognised “best practices” to establish the recommended testing requirements. Rigel Medical, their agents and distributors, accept no responsibility for any error or omissions within this booklet or for any misinterpretations by the user.

For clarification on any part of this booklet please contact Rigel Medical before operating any test instrument.

No part of this publication shall be deemed to form, or be part of any contract for training or equipment unless specifically referred to as an inclusion within such contract.

Rigel Medical assumes that the readers of this booklet are electronically and technically competent and therefore does not accept any liability arising from accidents or fatalities directly or indirectly from the tests described in this booklet.

Authors: Katherine Summers MEng and John Backes MA.

## 1 Introduction

It is estimated that 80% of hospitalised patients receive intravenous (IV) therapy and infusion devices are used extensively in clinical settings and patients' homes as an essential tool for providing perioperative care, critical care and pain management. The infusion of fluids uses a variety of designs providing the ability to feed, hydrate, medicate or replace blood loss [1-4].

An infusion system is the process by which an infusion device delivers fluids, nutrition and medication to a patient in a predetermined and consistent manner. They are capable of delivering medication such as insulin or

**Figure 1: Infusion device**



hormones, antibiotics, chemotherapy drugs, pain relief and can even be used for feeding (Fig. 1) [5].

An infusion pump is an electronic device used to control the administration of intravenous fluids to deliver measured amounts at careful and regulated rates. They often incorporate a mechanism or structure that mediate active transport across a biological membrane, using needles, where it has the best, most immediate effect. Occasionally subcutaneous, epidural or enteral methods are used but the amount of fluid is restricted, to as low as 3ml for subcutaneous methods [5, 6].

Infusion pumps have been widely used in clinical practice, such as in infusion IV therapy to infuse fluids, medication or nutrients into patient's circulatory system and can produce quite high but controlled pressure so as to inject controlled amount of fluids, however, pressure values vary with different pumps [7].


## 2 Introduction and History of Infusion

The prescription and administration of medication is an imprecise art as many dosages, for instance paracetamol, are not based on body weight and have only a maximum daily limit independent of size or weight of the user. On the other hand, treatment delivered through infusion directly into the bloodstream needs to be accurately controlled and often includes calculations based on the patient's physical condition (Fig. 2).

Infusion pumps can administer fluids in cases where other methods would be impractical or unreliable. For example, they can administer fluids at flow rates of 0.1 mL/hr (Alaris and CME™) which would be too small for a drip; provide medication at predetermined intervals (i.e. every minute); provide repeated boluses requested by the patient up to maximum per hour (i.e. patient-controlled analgesia); and administer fluids where volumes vary by the time of day and also for extended durations.

Although attempts at intravenous medicine date back to 1492, developments gained momentum in the 17th century with the first IV infusion device invented by Christopher Wren in 1658. However, developments were halted until the 19th century due to restrictions on IV infusion following deaths associated with blood transfusion after the first successful transfusion in 1665. During the 19th century the key elements of IV transfusion, which are still observed today, were established; a slow infusion process, awareness and prevention of risks from air embolism, and avoiding volume overload.

**Figure 2: Intravenous (IV) therapy**




One of the major developments in infusion pumps was the invention in the early 1970s of a wearable pump now known as the ambulatory pump by Dean Kamen, which enabled patient mobility during treatment [5].

There are many different types of infusion pumps, which are used for a variety of purposes and in a variety of environments; from simple gravity controllers which use a clamping action to vary the flow with the force of gravity, volumetric pumps which employ a linear peristaltic pumping mechanism and syringe pumps which work by pushing a plunger at a predetermined rate [4].

Occlusion is an obstruction or closure of a passageway or vessel. A blockage in infusion devices causes pressure to build up which can reduce the flow and cause harm to the patient. Most infusion pumps have a pre-determined occlusion pressure threshold See 5.3 [7], where an alarm activates once the pressure exceeds this limit. If the occlusion pressure alarm is set too high, the harmful effects can be prolonged prior to alarm. Therefore testing the occlusion pressure is crucial to infusion pump safety.

Infusion pumps must be periodically tested by qualified personnel to determine whether they are functioning properly [7]. There is a wide range of methods used to test the performance and accuracy of infusion devices which vary in procedure and equipment. However, the primary aim is to accurately measure the delivery volume and flow rate of the infusion device, check occlusion alarms and determine that it is safe for use.

**Figure 3: Infusion circuit**



An equally important point to testing infusion devices is not necessarily the accuracy of the Device Under Test (DUT) but how easy the test equipment is to operate, how simple it is to setup, if it can be left to run unattended and to ensure it produces the required information which is meaningful and easy to interpret. Therefore many users are moving away from basic measurements, such as weighing scales and burettes where continuous user input is required to ensure accuracy and moving towards automatic flow analyzers which record real time results graphically, such as the Multi-Flo from Rigel Medical.

### 2.1 Infusion Circuit

IVs are most often administered through bottles or bags of fluid that come premixed and range in standard sizes from 50 to 1000ml. The bag is hung from an IV pole and the IV tubing is attached to the bottom of the bag (Fig. 3). The IV tubing contains the drip chamber which is where measurements to calculate speed of a manual IV are taken. The drip chamber must always be half full to allow drops to be measured but also to prevent air from entering the tubing. The roller clamp is used to control the rate at which the IV fluid infuses. It can squeeze the tubing tightly, making it narrow and reduce the flow rate or it can loosen its pinch on the tubing to allow fluid to flow at a faster rate. Roller clamps should be closed before attaching the IV fluid to ensure no air gets into the tubing [6].

The safe range of daily fluid intake will vary based on the patient's condition, size and age. Several diseases such as heart failure, kidney failure, and

diabetes, for example, in a normal healthy adult, would require a total fluid intake (coming from IV fluids and/or fluids that they drink) of approximately 35-50ml/kg body weight/day. Also, IV administration is faster than any other method of administration because it goes directly into the blood, so it may be used when rapid action for treatment is necessary [6].

### 3 Types of Infusion

There is a variety of different types of infusion devices which all have individual procedures and mechanisms. However, they share an aim to provide accurate infusion at a predetermined rate over a set period of time. Infusion pumps in general use positive pumping actions which provide an accurate flow of fluids or drugs and all infusion pumps have commonalities including the alarm systems and the control panel etc.

Pumps are designed for a variety of clinical applications and their performance characteristics vary depending on the delivery volume, long & short-term accuracy and speed of the desired infusion. The Medicines and Healthcare Products Regulatory Agency (MHRA) have produced three categories in accordance with the potential infusion risks to aid selecting the most appropriate pumps for specific requirements (Table 1 overleaf).

The reliability of medical devices such as infusion pumps is extremely important because these devices are being used on patients who are likely to be in a critical condition [8]. Therefore, they incorporate warnings and alarms such as air in tube, excessive pressure, incorrect syringe size, loading problems, occlusion and also low battery warnings.

**Table 1: MHRA Categories of Infusion**

Therapy Category	Therapy Description	Patient Group	Critical Performance Parameters
A	Drugs with narrow therapeutic margin	Any	Good long-term accuracy Good short-term accuracy Rapid alarm after occlusion Small occlusion bolus Able to detect very small air embolus (volumetric pumps only) Small flow rate increments Good bolus accuracy Rapid start up time
	Drugs with short half-life	Any	
	Any infusion given to neonates	Neonates	
B	Drugs, other than those with a short half-life	Any except neonates	Good long-term accuracy Alarm after occlusion Small occlusion bolus Able to detect small embolus (volumetric pumps only) Small flow rate increments Bolus accuracy
	TPN Fluid maintenance Transfusions	Volume sensitive except neonates	
	Diamorphine	Any except neonates	
C	TPN Fluid maintenance Transfusions	Any except volume sensitive or neonates	Long-term accuracy Alarm after occlusion Small occlusion bolus Able to detect air embolus (volumetric pumps only) Incremental flow rates

There are five main types of infusion with the simplest device type being gravity controllers, which employ a clamping action to vary the flow of liquid due to the force of gravity. More complex systems use a positive pumping action for infusion; a simple version being an elastomeric pump which has a balloon reservoir that contracts to deliver fluid at a constant rate. Volumetric pumps may employ a linear peristaltic pumping mechanism applied to the infusion tubing or use a special cassette. Syringe pumps work by pushing the plunger of a disposable syringe along at a predetermined rate. The type of pump used is dependent on the patient's needs such as the required volume and the speed of the desired infusion [4].

### 3.1 Elastomeric Pumps


Elastomeric pumps, also called balloon pumps, are commonly used to administer liquid drugs such as local anesthetics or antibiotics. Fluid is held in a stretchable balloon reservoir and pressure from the elastic walls drive the fluid delivery with relatively consistent pressure until near the end of the infusion, which then spikes resulting in increased flow rate. The flow restrictor, usually a glass capillary or steel cannula, is molded into the tubing or situated within the reservoir to control the accuracy of the flow rate (Fig. 4).

Elastomeric pumps run independently without any electronics and are not gravity driven which makes them maintenance free and typically single-use. They are known for their reliability and accurate flow rate and drugs can be delivered




over a long period of time, up to seven days, making it ideal for outpatients or patients that require mobility. However, the pump's placement either above or below the patient will affect the flow rate.

**Figure 4: Elastomeric Pump**



**Figure 5: Syringe pump**



### 3.2 Gravity Controllers

For gravity controlled pumps the IV bag is situated higher than the patient's heart and fluid is gravity fed towards the patient using a roller clamp to control the flow. Timing how fast the drops fall through the tube drip chamber enables a flow rate to be determined and the greater the bag height, the quicker the flow rate in milliliters per drop (mL/drop). Each 300mm increase in height of the IV bag, produces approximately 25mmHg of pressure. This method of infusion is often used in Emergency Rooms, by paramedics in ambulances and by athletic trainers [4].

### 3.3 Syringe Drivers

Syringe drivers utilise an electronically controlled, electric motor to slowly depress the plastic syringe piston/plunger to drive the fluid into the patient at a predetermined rate. The electronics control the speed (flow rate), distance (volume infused) and the force (pressure) that the syringe plunger is pushed (Fig. 5). The operator must use the correct model and size of syringe, ensure it is correctly positioned and frequently monitor the delivery to maintain the expected drug dosage up to 100ml of fluid at flow rates between 0.1 to 100mL/hr.

Syringe drivers are the preferred choice for lower volumes and low flow rates and can infuse small volumes of fluid over an extended period of time while maintaining a constant rate. However, users should be aware that the flow delivered at the start of any infusion might be less than the predetermined value due to mechanical slack which must be taken up before a steady flow rate is achieved. At low flows it can be some time before any fluid is delivered to the patient.

### 3.4 Volumetric Pump


The volumetric infusion pump is the most sophisticated type of pump which forces fluid into the patient's vein under pressure and resistance [6]. The IV bag is situated higher than the pump and patient and the pump utilises a linear peristaltic action or uses a special cassette, known as a piston cassette, to control the infusion fluid. A set of rollers pinch the flexible tubing to push the fluid towards the patient (Fig. 6). Volumetric pumps are used to administer up to 1000mL of fluid at flow rates between 0.1 to 1000mL/hr. However, this type of pump is not considered accurately appropriate for the delivery of fluids at a rate lower than 5mL/hr.

Peristaltic action is a continuous rippling wave motion which can be linear or rotary. The Baxter Colleague™ pumps use a linear peristaltic pumping action and the rotary peristaltic action is most commonly used in feeding and Patient Controlled Analgesic (PCA) pumps.

Cassette type pumps have a plunger mechanism built into the tubing set where the stepper motor drives the motion. The cassette also has a valve which is timed to only open when the plunger is

pumping. This valve can sometimes be used as a free-flow mechanism. The Ivion Kids pumps used this type of pumping mechanism as do some PCA pumps.


**Figure 6: Volumetric pump**



### 3.5 Patient Controlled Analgesia (PCA) Pump

A PCA pump is a type of syringe driver that allows the patient, within specified limits, to control their own drug delivery through a patient hand control, which when pressed, delivers a pre-set amount of analgesic drug, known as a Bolus (Fig. 7). After the delivery of the bolus the pump will return to the normal basal flow rate and refuse to deliver another bolus until a defined time has passed, known as the lockout time, thus limiting the number of boluses allowed. The bolus size and lockout time, along with the constant basal drug infusion rate are pre-programmed by the clinician.

**Figure 7: Patient controlled analgesia (PCA) pump**




### 3.6 Ambulatory Pumps

Ambulatory pumps are small, light, battery powered syringe or cassette mechanisms which are designed to be portable or wearable giving the patient freedom to move when receiving treatment. This allows for medication to be administered on an outpatient basis, especially for those who need round the clock injections (Fig. 8) [5].

Most ambulatory pumps have minimum alarms, which mean that patients and carers need to closely monitor the infusion. Consideration also has to be given for the hazards a portable device is exposed to e.g. knocks, fluids, electromagnetic interference etc. In general, critical drugs which require a constant flow should not be administered using ambulatory pumps.

**Figure 8: Ambulatory pump**



### 4 Infusion Characteristics


In an infusion system, the pump uses pressure to overcome the resistance to flow to deliver infusion. The greater the resistance in the IV circuit the higher the pressure is required to deliver the prescribed flow. Resistance arises due to filters, anti-siphon and anti-reflux valves, administration set, in particular the cannula, the internal diameter and potential kinking of tubing. This along with viscous solutions and syringe/cassette stiction can accumulate the pressure [7]. The pressure available is related to the height, that the bag or bottle of fluid is situated above the patient's heart, with increased height produces increased pressure. Therefore infusion devices must be capable of delivering infusions at pressures between 100 and 750mmHg (2 to 15PSI) to overcome any internal or external resistance to the flow of fluid.

The patient's venous system contains pressure and therefore the IV infusion device must be proven to be compatible with pressures found in the body. To deliver simple electrolytes into an adult intravenous site at 100mL/hr the pressures associated with resistance are summarised below (Table 2).

**Table 2: Resistances within the body to infusion**

Point of Resistance	Pressure (mmHg)
Maximum adult venous pressure	30
Filter	10
Cannula	100
Administration set	1
Total = 141	

**Figure 9: Fluid delivery pattern for volumetric, syringe and ambulatory pumps**



The pattern of fluid delivery is dependent on the type of pump used. The typical flow patterns for volumetric, syringe and ambulatory pumps at a flow rate of 1mL/hr (Fig. 9). Each pump delivers fluid accurately at 5% of the set flow rate over long periods. However, only the syringe pump delivers fluid accurately over very short periods of time.

### 5 Problems Associated with Infusion

Millions of infusion devices are used in hospitals and in the community every year. The vast majority of infusions are delivered safely. However, incidents with infusion devices continue to dominate adverse incident reports with at least 1000 incidents investigated by the

MHRA between 2005 and 2010 in the UK alone [4, 10]. In the past five years, the FDA has received more than 56,000 reports of adverse events associated with the use of infusion pumps. Those events have included serious injuries and more than 500 deaths. Between 2005 and 2009, 87 infusion pump recalls were conducted to address identified safety concerns, according to FDA data [14]. The majority of problems relate to over-infusion of drugs, either due to user error with dosage and patient data, product design and engineering or software malfunction. In many fatal incidents no fault has been determined with the infusion device which suggests that user error is a contributing factor whereas, in practice it seems that syringe pumps produce the most significant problems compared to other methods of delivery [4, 8].

Genuine infusion pump malfunctions rarely occur, except when a pump has been mishandled, dropped or damaged. However, faulty or mismatched batches of administration sets can cause problems. Prompt reporting of damaged pumps, suspect administration sets or unexplainable alarms can lead to safety checks and investigations which improve the overall drug infusion safety [8].

There are many dangers when using infusion devices which need to be minimised for patient safety and therefore most infusion devices incorporate detecting warnings and alarm systems including: air in tube, excessive upstream pressure and both up and downstream occlusions, syringe empty / nearly empty and also low battery signals.

A careful set-up procedure will prevent the majority of errors; however, major problems can still occur after the onset of an infusion. This is particularly the case when patient care is transferred to different staff or clinical area. A proper hand-over with complete documentation is necessary. Good communication such as a checklist, is the most important factor in reducing infusion risks and increasing compliance with correct procedures [9].

### 5.1 Air Embolism

The problems of air entrapment are restricted to volumetric pumps and therefore IV pumps are recommended to have a means of detecting and preventing air in the tubing using either a photo detector or ultrasonic sensor which is capable of detecting single bubbles of about 0.1mL. Many systems of administration employ a drip chamber, which prevents air from entering the blood stream, known as an air embolism [9].

The ultrasonic air detection system responds to any air in its pathway and if tubing is not correctly attached it may give a false alarm. Also, solutions which are prone to froth can activate the sensor.

False alarms are at best a nuisance as air in the tubing can be dangerous to the patient if it gets into the blood as approximately 3-8mL of air per kg can cause cardiac arrest, that's about 210-560mL for a 70kg patient. Even if there appears to be no air in the line, the action of pumping can often draw air out of the solution and this can accumulate into significant sized bubbles or if there is a leak in the upstream line, air can be drawn into the line. Air entrapment can be prevented by the intrinsic design of the pump or by using an air detector.

### 5.2 Free Flow or Siphonage

If the pump is higher than 300mm above the patient's heart, and the roller clamp is fully open on a basic gravity infusion set-up this will allow free-flow or siphonage which allows all the fluid and any air in the infusion bag to infuse into the patient. Siphonage should not occur if the administration set for the appropriate infusion device is loaded correctly and the following simple precautions have been taken:

- The tubing, syringe or cassettes are correctly loaded.
- The pump door is closed.
- Syringe pumps should be placed at the level of the patient with the syringe securely located in the mechanism and with an anti-siphon valve in the line.
- The absence of free-flow can also be checked by opening the roller-clamp before connecting to the patient.

Confirmation that the roller clamp is secure before removing the administration set from the pump is essential as this is a primary means of occluding (causing a blockage) in the line. If the administration set is fitted with a flow-stop mechanism, this should not be solely relied upon and the roller clamp should still be used [7, 9].

### 5.3 Occlusion

An infusion pump attempts to maintain sufficient pressure on the fluid to cause it to flow at the set rate and detects resistance to flow, allowing the pump to increase the pumping pressure to maintain the set flow rate. Any occlusion or blockage can result in patient harm caused by

increased pressure in the line and an interruption to the therapy. This is typically due to a partial or complete block in the delivery tubing e.g. kinks in the tube, clamp or tab closed; or in the cannula e.g. clotted off or a change in position. Occlusion can cause an interruption to the delivery and effectiveness of the therapy. Therefore, occlusion alarms are used to indicate when the pump is unable to sustain the set flow rate or when pressure in the line increases, which is usually a stall sensor on the IV pump stepper motor or a pressure transducer within the line. Any blockage can cause the pressure in the line to increase which causes the tubing to expand and if the pump detects an unacceptably high pressure, known as the occlusion limit, the pump's occlusion alarm will sound. Adult alarm settings should be approximately 150mmHg above working pressure and therefore a default setting of 300mmHg is standard [7].

#### 5.4 Bolus

A syringe that is placed into a pump whilst connected to the patient is highly likely to deliver an infusion bolus. Syringes should always be inserted into pump mechanisms before being attached to patients and the infusion line should be temporarily disconnected. Even if the infusion set is left connected to the syringe driver or volumetric pump during transfers, large alterations in the device height relative to the patient may lead to over or under-infusion because of tubing compliance [7].

Once occlusion is cleared a sudden bolus will infuse into the patient, which is the additional fluid generated in the line during the occlusion.

This will cause a momentary raise in the fluid's concentration. The higher the occlusion alarm/pressure limit is set, the longer IV related problems are allowed to develop before the alarm alerts staff of any issue and resulting in a larger bolus. Flow rate also affects alarm response time as higher flow rates result in more rapid alarm response than lower flow rates.

To prevent this bolus occurring, the pressure in the system should be reduced by temporarily opening the system to air or disengaging the clamp on the syringe plunger. Modern syringe drivers have an automatic 'back-off' facility when occlusion is detected which briefly draws back the plunger thus reducing the size of any bolus [7, 9].

#### 5.5 Vein Closing After Infusion

All IV pumps have a 'keep vein open' (KVO) function which occurs when the time or volume to be infused is reached and prevents the infusion needle from clogging by continuing to infuse a flow between 1 and 5mL/hr, which is always a rate lower than the set infusion rate [7, 9].

#### 5.6 Infiltration or Tissueing

Infiltration or tissueing is the leakage of fluids as a result of improper infusion and can be caused by punctures of vein wall, dislodgement of catheter/needle cannula, inappropriate venous site, improper cannula size or excessive delivery rate.

Infiltrations can occur at pressure greater than 100mmHg and pumps deliver infusions at pressures of 750mmHg (15PSI) to overcome fluid resistance which makes them capable of

producing extra-vascular patient damage. To minimise this risk infusion pumps have adjustable pumping/occlusion pressure settings (normally technical settings), newer units incorporate user settings.

A balance is required between setting the pumping pressure to overcome all baseline and intermittent resistance whilst preventing the possibility of tissueing. Setting the pressure low, at best causes nuisance occlusion alarms.

## 6 Testing Infusion Devices

### 6.1 Introduction

Optimal infusion is the ability of a device to reliably deliver the prescribed dosage and volume to the patient, at a pressure which overcomes all baseline and intermittent resistance, whilst causing no harm to the patient [10]. The reliability of infusion pumps is extremely important because these devices are used for patients who could be in a critical condition; also due to incidents associated with infusion devices there is a need to adequately validate the accuracy and performance of such devices [7].

All tests carried out should mimic real life settings, as close as possible, which are carried out by nursing staff when the system was initially set up. Hence, testing should reflect what the manufacturer recommends ensuring that the equipment is working within its specification. The user should always use the specification for reference when taking measurements and/or use the test equipment recommended or specified by the manufacturer.

There is a variety of methods currently used to test infusion devices and determine their performance accuracy. They vary in procedure and equipment but fundamentally, the aim is to measure the accuracy of the delivery volume and flow rate over a range of time periods, typically between 10 minutes and 1 hour.

Common flow measuring principles:

1. Volumetric - Flow calculated after a certain volume has been delivered. The greater the volume over a certain time, the greater the flow.
2. Mass - Flow is calculated based on temperature difference between two points within the sensor, the greater the temperature difference the lower the flow.
3. Bubble tracking - Flow calculated based on the displacement of an inserted airbubble into the flow sensor part. The greater the displacement, the greater the flow.
4. Pressure based - Flow is regulated within the flow sensor to a set line pressure. The greater the potential pressure built-up in the line, the greater the flow rate.
5. Displacement of syringe plunger - Flow rate is calculated based on volume displaced by the syringe plunger over time. The syringe type and volume are required to provide an accurate calculation.

Occlusion and alarm pressures and also bolus delivery must also be tested to maintain the

performance of the infusion device especially in PCA devices where the bolus is self-medicated. A visual inspection and electrical safety test should also be considered to make sure all aspects of patient safety and instrument reliability are covered during the test procedure.

When looking at the accuracy of an infusion system the system must be looked at as a whole, and not just the pump or driver individually as the inaccuracies within the measuring device need to be added to the accuracy limits of the DUT. Even though the infusion device itself may only have a small percentage error; taking into account all possible inaccuracies of the syringe used and other external equipment, including the tubing set, could increase the inaccuracy to 10%. For example, the diameter, length and elasticity of the tubing can provide greater inaccuracies than the pump itself, especially over time as the material's properties change. Due to these inaccuracies it is generally accepted that the testing method and equipment used must be an order of magnitude of 1/10th more accurate than the device being tested [8].

## 6.2 Volume/Flow Rate

The most commonly used methods for measuring volume or flow rate are direct volumetric; using either graduated cylinders or burettes; and, derived mass measurement using a measurement vessel and weighing scales. When using these methods it is recommended that only distilled water be used to measure the total volume delivered as the weight in grams is equivalent to the weight in milliliters.

### 6.2.1 Cylinders and Burette

When using a direct volumetric method the estimated time is calculated prior to testing by first determining the rate per minute and then the desired volume over the rate per minute, to determine the time to reach the volume.

For example:  $q=V/t$

When testing a pump with a flow rate of 240mL/hr, the rate per minute =  $240/60= 4\text{mL}$  per minute. If the desired measurement volume is 20mL, the desired volume/ rate per minute =  $20/4 = 5$  minutes. So, it would take 5 minutes to reach 20mL of volume with a flow rate of 240mL/hr.

The pump is set to run at the desired flow rate and the amount delivered into the measuring device is measured and the total time recorded, often using a stopwatch. The results are then used to determine the accuracy of the system and the average rate can be calculated as a total volume over total time (mL/hr).


Some technicians produce test rigs which include a connection between the stopwatch and the pump so that when the pump finishes delivering fluid the stopwatch automatically stops. Otherwise the user must continuously monitor the test to record the volume at the desired time.

Reading the volume from a measuring cylinder accurately is difficult as the eye must be in line with the lower level of the surface of the water, known as the meniscus, and the measuring cylinder must be placed on a level surface [12].



Also when using a measuring cylinder the unit graduations must be considered which provides a degree of inaccuracies e.g. a 100mL cylinder with 1mL graduation will therefore only be accurate to  $\pm 0.5\text{mL}$ . The total capacity of the cylinder must be considered not just the volume filled (Fig. 10a).

**Figure 10(a):**  
**Measuring cylinder**      **10(b): burette**




Burettes are calibrated, marked in divisions, and are usually more accurate than measuring cylinders. The amount of liquid that is tapped from a burette can simply be read from the burette by checking the level of the meniscus

(Fig. 10b). A burette is particularly suitable for measuring volumes less than 50mL with a tolerance of  $0.1\text{cm}^3$  [14]. Uncertainty in the measurement of a 50mL burette is 0.1mL therefore  $\pm 0.2\%$ . Adding this error to the typical acceptable tolerances for a syringe driver of 2% will give the total inaccuracy of  $\pm 2.2\%$ . The total inaccuracy needs to be checked against that of the manufacturer's recommended limits to make sure that it is within allowable tolerances.

### 6.2.2 Weighing Scales

One of the most accurate and simplest methods to measure volume and flow rate is to use scales to determine the weight or mass of the infused liquid (Fig. 11). Scales can be calibrated to provide extremely accurate results and are useful for measuring volumes as the weight of the empty cylinder can be "zeroed out" (tared) hence the results are a direct reading of the weight of the fluid volume delivered based on the assumption that 1g is equivalent to 1mL of fluid (derived mass method).

**Figure 11: Weighing scales**



The pump is set to deliver a specific volume at a specified rate and the vessel is re-weighed at the end of the test. A stopwatch would be required to record the time. This method requires user input to monitor the time and to maintain accuracy. A percentage error can be calculated from the figures obtained.

### 6.2.3 Vernier Calipers/Dial Gauges


Vernier calipers and dial gauges give a direct reading of the distance measured with high accuracy and precision. They are functionally identical, with varying ways of reading the result. Calipers are comprised of two jaws, one attached to a fixed scale and the other attached to a sliding scale. In order to measure an object's width, the object is simply placed between the caliper's jaws and then the distance between the jaws is read to an accuracy of 0.05mm (Fig. 12a) [11].

Dial gauges are a precision measurement capable of producing extremely fine measurement increments up to 0.001mm. The measurement inputs are transferred to the gauge via a plunger or hinged lever and are available with either an analogue needle and dial indicators or digital liquid crystal displays (LCDs) (Fig. 12b) [10].

These methods are used to check the plungers travel accuracy. Some pumps such as the Alaris (Carefusion) pumps come with a dedicated test set that uses a linear test gear and a stop watch to determine how far the plunger moves to check incremental and overall time of delivery.

**Figure 12(a):  
Vernier calipers**

**12(b): dial gauge**



### 6.2.4 Electronic Devices

Automated analyzers allow the user to set up a test and the test will run unassisted thus allowing the user to move away from the device during testing. Electronic analyzers provide real-time records of the delivery rate and volume allowing for continuous infusion device testing without constant supervision. Analyzers provide continuous spot readings which can then be graphed to provide quick to read results; where problems and anomalies are easily visible. Electronic calibrated analyzers provide automation for batch and multi-channel testing with some analyzers having up to 4 independent channels which can run simultaneously, reducing overall test time while maintaining a good degree of accuracy and minimising user input and any errors associated.

These devices work by timing how long a small internal volume takes to fill and therefore measures the volumetric accuracy of the infusion device where accuracy is determined by the measuring chamber accuracy and its resolution is determined by the chamber volume. It provides the user with a method to obtain a quick check as to whether a pump is generally working. You can test it over a range of flow-rates and time periods with consistent accuracy [9]

### 6.3 Occlusion and Alarm Pressure

Occlusion is the interruption or ineffective therapy of infusion due to a blockage, momentary closure or obstruction of the passageway. Occlusion sensors can detect both upstream and downstream restrictions and the downstream occlusion sensor sensitivity can be adjusted to suit the needs of the patient/hospital through the pump configurations.

Dynamic monitoring systems provide the ability to monitor downstream pressure or resistance allowing for rapid detection of full or partial blockages. To prevent nuisance alarms and interruptions to therapy the initial occlusion pressure needs to be set above the systems normal running pressure.


Occlusion alarms are used to indicate when the pump is unable to sustain the set flow rate and therefore pressure in the line begins to increase. This is typically due to a partial or complete block in the delivery tubing e.g. kinks in the tube, clamp or tab closed; or blockage signs in the cannula e.g. clotted off or a change in position.

#### 6.3.1 Analog Pressure/Force Gauge

Gauge pressure sensors measure the pressure relative to the atmospheric pressure which surrounds it. The term pressure gauge usually refers to a self-contained indicator that used flexible elements as sensors which converted the detected pressure change into the mechanical motion of the flexible element which then rotates a pointer in front of a dial, see (Fig. 13). In these early mechanical pressure sensors, a Bourdon tube, a diaphragm, or a bellows element detected the process pressure and caused a corresponding movement.

Automatic control systems have since evolved and the free end of a Bourdon tube (bellows or diaphragm) no longer has to be connected to a local pointer, but served to convert a process pressure into a transmitted (electrical or pneumatic) signal. At first, the mechanical linkage was connected to a pneumatic pressure transmitter, which usually generated a 3-15psig output signal for transmission over distances of several hundred feet, or even farther with booster repeaters. Later, as solid state electronics matured and transmission distances increased, pressure transmitters became electronic. The early designs generated DC voltage outputs (10-50mV; 1-5V; 0-100mV), but later were standardized as 4-20mA DC current output signals. [12].

**Figure 13: Force dial gauge**




A force gauge can be used where the occlusion pressure is measured in kgF and for a set occlusion pressure 'level', the manufacturer should state that the pump should alarm when a certain amount of force in KgF is applied to the plunger, see (Fig. 13), with stated limits.

### 6.3.2 Digital Pressure Meter

Digital pressure meters (Fig. 14) are accurate and can be used to carry out the occlusion testing

where the pressure is increased by the pressure meter until the alarm is activated. The pressure of the meter at which the alarm sounded is recorded and compared to the alarm pressure set on the infusion device to determine accuracy. This method produces good, accurate and repeatable results, however it lacks the ability to produce a report so all results are recorded manually.



**Figure 14: Digital pressure meter**

### 6.3.3 Electronic/Automatic Devices

Electronic devices will use a pressure meter to increase the pressure in the infusion device which is intended to pass the occlusion limit and set off the occlusion alarm.

With manual occlusion pumps, often the user must indicate on the analyzer when the alarm is activated to record the occlusion pressure at alarm. Some modern infusion devices have an auto rewind function where the pump mechanism reverses when occlusion is observed. This is to try and reduce any ill effects caused by the increased pressure in the line. Automatic analyzers will display the maximum occlusion pressure at which auto rewind occurs without the need for user input.

### 6.4 PCA

A PCA test is to check the bolus which can be administered by the patient and to make sure that the correct volume and also lockout time are working so that the patient cannot over infuse medication to themselves. To test the lockout facilities on PCA pumps, a time is set; e.g. 5 minutes and then a bolus is demanded and then subsequent demands are made periodically to establish when the next bolus is permitted. This should be then repeated at a variety of volumes such as 0.1, 1.0 and 5.0mL target bolus volumes to maintain consistency [12].

#### 6.4.1 Electronic/Automatic Analyzers


With electronic infusion analyzers the user will set up the initial basal flow rate and the bolus volumes to be delivered. The infusion will begin at the set flow rate and when bolus is activated the

analyzer will measure the change in flow rate and record the bolus over time, often producing a real-time graphical display of the increased volume and flow rate until the desired bolus is infused and the flow returns to the basal rate. Boluses can then be repeated to make sure the lockout time on the infusion device does not allow for repeat boluses within a specific time. Often five boluses are delivered and the average taken to determine the performance accuracy of the pump.

### 6.5 Trumpet Curve

Trumpet curves show the accuracy performance of an infusion device at set intervals in the second hour of infusion and during the final hour of infusion as required by IEC 60601-2-24. The trumpet curve indicates the maximum percentage deviation, both positive and negative, from the expected flow rate relative to the time interval, known as the observation window. For example; for one hour the overall deviation may be -2%, whereas over an interval of two minutes the deviations can vary between +7 and -10% (Fig. 15) [13].

Initially the pump needs a period of time to 'warm up', this is often referred to as the settling in period where the flow rate can vary as the infusion settles. The flow rate of infusion devices is defined as settled after one hour of infusion and the majority of pumps have a stated accuracy of  $\pm 3\%$  of the set flow rate after this one hour settling period. Therefore, trumpet curves are produced during the second and occasionally the final hour of testing to determine whether the manufactures specified performance accuracy is met [13].



**Figure 15: Example of a trumpet curve**

Some manufacturers have begun to include trumpet curves in their user manuals. The graph has developed its name as it looks like a trumpet which converges to the right side. Time units are positioned on the x axis and percentage on the y axis (Fig. 15). These  $E_p$  (max) and  $E_p$  (min) curves represent the maximum positive and negative percentage deviation from the expected dose for a given observation window during the infusion.

As the accuracy of a medical device also depends on the set flow rate, often a range of curves are included. The user can then correlate the half-life of the administered fluid with the observation interval to determine the suitability of the pump for a particular application.

However, in general terms, the longer the time interval, the more accurate the dose [13].

## 7 Multi-Flo

### 7.1 Introduction

The market defining Rigel Multi-Flo Infusion Pump Analyzer is a portable instrument to accurately and swiftly verify the performance of all infusion devices.

Accurate and fast measurement of flow (from 10 $\mu$ L to 1500 mL/hr), pressure, PCA (bolus) and total volume ensure all requirements are tested according to IEC 60601-2-24 during the manufacturing, design and performance test stages.

On-board memory stores test data and allows fast transfer to the PC for traceability. Fast and error-free asset information can be directly entered into the Multi-Flo via a compact Bluetooth barcode scanner or an optional USB keyboard.

An upgraded and enhanced version of Rigel's Med-eBase software provides the ability to create test templates, custom test certificates and the ability to control and configure the Multi-Flo infusion pump analyzer from a computer.

Remote control allows the user to create definable tests and test sequences. The Multi-Flo can then run automatic tests from the built in sequences or create custom and/or manufacturer specific test sequences using Med-eBase software.

**Figure 16: Four Channel Multi-Flo Infusion Analyzer**




PC control provides real time results on the PC screen and can be stored against the device under test (asset number) in the database. The results can be viewed in real-time and after the test is completed.


## 7.2 Testing Infusion Devices with the Multi-Flo

The Rigel Multi-Flo is easily connected to both syringe driver and volumetric pumps as shown in the diagrams below. Ensure the flow direction is as per diagrams. The flow inlet is always the top connection and the flow outlet is positioned below the inlet for each channel.

### 7.2.1 Syringe Driver Connection



### 7.2.2 Volumetric Pump Connection





### 7.2.3 Volume/Flow

The Rigel Multi-Flo is capable of measuring instantaneous flow at a resolution of 10µL/hr. In addition, the flow rate can be viewed based over an average, user selectable, period. Detecting peak and minimal flow rates on real time curves, the flow measurement provides the benefit of

faster test times at low flow rates. The ability to detect low flow rates makes the Multi-Flo infusion pump analyzer a versatile tool for all types of infusion. Custom tests and sequences can be created on remote control software and user-definable limits help clearly indicate whether the performance is within the manufactures specification.

**Figure 17: Volume/Flow rate test setup, results screen and real-time graph from Med-eBase software**






### 7.2.4 Occlusion




The Multi-Flo Occlusion test simulates an obstruction in the infusion process and monitors the variation in pressure due to the blockage. Most infusion devices have the ability to detect this obstruction and provide an occlusion alarm. The occlusion test is able to test this alarm feature in infusion devices.

Some infusion devices have an auto rewind function at the occlusion alarm and the Multi-Flo can detect the maximum occlusion pressure at which auto rewind occurs. In infusion devices that do not have auto rewind (normal stop/alarm pumps) the user must press the Green button on the Multi-Flo when the occlusion alarm sounds which will stop the test then record and display the maximum occlusion pressure observed.

**Figure 18: Remote control occlusion set up and test screen**



**Figure 19: PCA set up and test screen**



### 7.2.5 PCA

The PCA test determines the additional volume delivered on top of the basal flow rate set by the user. The additional volume, sometimes referred to as bolus, is an indication of the correct safety settings of an infusion device. The user needs to enter the basal flow as the basal flow rate setting is used to determine the additional volume being delivered i.e. the bolus.

## 8 Conclusion

Whatever method chosen, whether it is with a burette (direct volumetric measurement), weighing scales (derived mass measurement) or automatic analyzers, the most important consideration is that the method provides reliable and accurate results.

Direct volume measurements are techniques which provide a good degree of accuracy without the need to worry about the design or performance accuracy of another piece of equipment. However, they are very labor intensive and require continuous user input.

Electronic, automated devices combine the various manual test methods with the ease-of-use of one analyzer that offers the ability to conduct the most important aspects of infusion device performance accuracy analysis with little user input and the ability to test multiple infusion devices simultaneously.

The volume of fluid delivered is used to measure the system's volumetric accuracy and other tests determine the performance accuracy of the device to maintain safety for staff and patients who will use the equipment. Most manufacturers specify system

accuracy under stated test conditions including the infusion set, temperature of test, rates, etc. this indicates that there are a number of elements that can affect the accuracy and therefore the device under test must be used within the recommended specification to maintain the given limits of accuracy.

Most adverse incidents are eventually identified as user error. Designing for usability can contribute to minimizing user error. Work has been initiated to develop a formal ergonomic testing procedure for application to all devices. At present, user instructions are assessed for clarity and readability, conciseness, and indexing. Procedures for using the device are systematically worked through on the bench after other testing is completed. Any hazardous potential misuse is noted. For all available alarms, the reliability, readability of text displayed, the alarm tone quality, the positioning of alarm lights, and methods of silencing alarms are tested as part of the ergonomic assessment of the device [12].

The technological advances in infusion pumps during the past forty years have transformed the treatment of patients in hospitals, as well as afforded the ability to infuse fluids on an outpatient basis or in a home environment and enable patients to receive treatment while going about our daily lives [5].

Therefore, measuring performance accuracy needs to keep pace with the changing and advancing technology of new pumps, which can infuse very low volumes and for extended periods, to be able to evaluate them in terms of volume/flow rate, occlusion and bolus measurements.

## 9 References

1. NHS, National Patient Safety Agency. (2004). *Improving infusion device safety. Safer Practice Notice*. Issue 1; 20 May 2004.
2. Meier, B. (2010). *F.D.A. Steps Up Oversight of Infusion Pumps*. Available: [http://www.nytimes.com/2010/04/24/business/24pump.html?\\_r=2&](http://www.nytimes.com/2010/04/24/business/24pump.html?_r=2&). Last accessed 27th August 2013.
3. Rosenthal, D. (2006). *History of IV therapy*. Available: <http://workplacenurses.com/id39.html>. Last accessed 25th August 2013.
4. MHRA. (2010). Infusion Systems Device Bulletin. *Safeguarding public health*.
5. David Ostendarp. (2010). *A Brief History of Infusion Pumps*. Available: <http://www.articleonlinedirectory.com/322515/a-brief-history-of-infusion-pumps.html> . Last accessed 25th August 2013.
6. Wladis, C. (2012). *Lecture 6: Introduction to IVs*. Available: <http://www.cwladis.com/math104/lecture6.php>. Last accessed 28th August 2013.
7. Zhang, P et al. (2009). Design of occlusion pressure testing system for infusion pump. *Biomedical Science and Engineering*. 2 (6), 431-434.
8. Davis, W. (2010). *Infusion Device*. Available: [http://www.ebme.co.uk/arts/infusion/inf\\_train.php](http://www.ebme.co.uk/arts/infusion/inf_train.php). Last accessed 27th August 2013.
9. Keay, S and Challander, C. (2004). The Safe Use of Infusion Devices. *Anaesthesia, Critical Care & Pain*. 4 (3), 81-85.
10. Paul Scott. (2013). *What Is a Dial Gauge?* Available: <http://www.wisegeek.com/what-is-a-dial-gauge.htm>. Last accessed 25th August 2013.
11. *How to use Vernier Calipers*. Available: <http://www.physics.ncsu.edu/courses/pylabs/caliperHowTo.pdf> Last accessed 25th August 2013.
12. Bath Institute of Medical Engineering. (2005). *Test Methods and Protocols*. Available: <http://www.bath.ac.uk/bime/evalcentre/publications/protocols.pdf>. Last accessed 28th August 2013.
13. Ferrari, R and Beech, D. R. (1995). 1. Infusion Pumps: Guidelines and Pitfalls. *Australian Prescriber*. 18, 49-51.
14. U.S Food and Drug Administration. (2010) FDA Launches Initiative to Reduce Infusion Pump Risks Available: <http://www.fda.gov/NewsEvents/Newsroom/PressAnnouncements/ucm209042.htm> Last accessed 19th May 2014.

Rigel Medical, Seaward Group USA,  
6304 Benjamin Road, Suite 506,  
Tampa, Florida, 33634, United States

Part of  **SEAWARD**  
GROUP USA

**Tel:** 813-886-2775 **Email:** [enquiry@rigelmedical.com](mailto:enquiry@rigelmedical.com)

**Web:** [rigelmedical.com/usa](http://rigelmedical.com/usa)

