



RIGEL MEDICAL

EYOYO Bluetooth® Scanner Pairing Guide

When pairing an EYOYO Bluetooth® scanner to your instrument, always save it in the dedicated Barcode field within the Bluetooth® pairing menu. Saving the scanner in any other pairing field may prevent a successful connection.

Ensure that in the *Menu* →  > *Setup* > *User Admin* > *Preferences*, the data entry field is set to Barcode. Failure to enable this setting will result in text entry only in the Asset ID field.

How to Pair the EYOYO Bluetooth® Scanner

- On your Rigel instrument navigate to the Bluetooth® favorites menu.
- Select the Barcode field then access the setup menu .
- Enable pairing mode on the scanner by scanning the two barcodes below. Scan the left barcode first, then scan the second barcode. The scanner will emit a beep after each successful scan, and the blue indicator light will flash when it is ready to pair.



- On the Rigel instrument, press the search button  to search for all available Bluetooth® devices.
- Once the search has completed, the scanner name is displayed as EY-XXXX-XXXX.
- Use the arrow keys to select this device from the list, then press the **+** button to add it to the Barcode field on your instrument.
- Press the **✓** button to confirm your selection.
- Press the **✓** button in the Bluetooth® favorites screen to save your settings.
- From the home screen, press the menu button →  then select the auto mode option and confirm your selection **✓**.
- Ensure that the Asset ID field is selected and that the scanner is powered on.
- The instrument will search for the scanner as seen by the hourglass indicator  on the top right corner on the screen.
- After a few seconds, the scanner will emit a beep and pair with the instrument. On the top right-hand corner of your screen a Bluetooth®  symbol will appear to indicate a successful connection.
- It is possible to now scan a barcode to add your Asset ID.

UK ☎ : +44 (0) 191 586 8730 | ✉ sales@seaward.com

USA ☎ : +1 732 287 3680 | ✉ sales@seaward-groupusa.com

Germany ☎ : +49 911 8602-0 | ✉ sales@gossenmetrawatt.com

rigelmedical.com