

Maximum Data on the Rigel 288

Max number of asset records	5,000
Max number of test sequences	50 (of which 12 are pre-defined)
Characters in a test sequence name	25
Characters in a User Test name	200
Characters in a User Test unit	25
Max number of user names held in list	40
Max number of entries in the comments string list	80
Max number of entries in the AP module name list	20
Number of Trace variables (Site, Location etc.)	8
Max number of entries in each trace variables list	40
Max number of tests per test sequence	100
End.	Max of 300 Test Codes depending on size of each test code. (If each code references lots of Test Sequences, then this number will reduce)