



Innovating Together

Rigel Uni-Therm Accessory List

Rigel Uni-Therm

Accessory List

Foot Switch Adapter Cables				
Part No.	Description	Compatibility	Image	Specification
367A950	Oscilloscope connection lead	All Manufacturers		Oscilloscope SMB to BNC connection lead
367A951	Foot switch adaptor	ERBE ICC350		4 pin male
367A952	Foot switch adaptor	Valleylab™		4 pin male
367A953	Foot switch adaptor	BOWA ARC 200, 300 and 400		7 pin male
367A955	Foot switch adaptor	ForceTriad™		6 pin male

Part No.	Description	Compatibility	Image	Specification
367A956	Foot switch adaptor	Conmed		4 pin male
367A957	Bipolar foot switch adaptor	Valleylab™, Covidien, Megadyne		3 pin male
367A958	Foot switch adaptor	Olympus PSD/JES 40 and 20		5 pin male
367A959	Bipolar foot switch adaptor	Olympus		3 pin male

Neutral Plate Lead Cables

Part No.	Description	Compatibility	Image	Specification
800A990	Split neutral plate lead with CQM	Aksa, Martin, Eschmann, Bovie (Aaron), Integra, Berchtold, Erbe, Aesculap		2 pin female with plastic REM control (0.3m length)
800A991	Single neutral plate lead 6.3 Jack	Valleylab™, Ellman, Integra Intern, Berchtold Intern		6.3 Jack (0.3m length)
800A992	Split neutral plate lead with CQM	Valleylab™		2 pin female and REM control (0.3m length)
800A996	Single neutral plate lead w/o CQM	Valleylab™		2 pin female (0.3m length)

Accessory Cables

Part No.	Description	Compatibility	Image	Specification
800A993	Foot switch cables	All manufacturers		4 x 1.5 metre cables 2 x 4mm plugs 2 x blue and 2 x yellow
800A994	Active cable	All manufacturers		1 metre red cable 2 x 4mm plugs
800A995	Bipolar cable	All manufacturers		1 metre black cable 2 x 4mm plugs
800A997	Accessory cable set	All manufacturers		Consists of: 800A996, 800A995, 800A994, 800A993, 800A992.
286A954	Crocodile clip	All manufacturers		Additional accessory which may be required during performance testing
367A954	Jumper links (per 3)	Rigel Uni-Therm		3 x 4mm jumper links for side panel of Rigel Uni-Therm

Data Communication

Part No.	Description	Image	Specification
396A976	USB download cable		USB cable for data communication with PC and performance analyzers

Software

Part No.	Description
383A910	Med-eBase V2 Allows you to manage your assets the way you want to
383A914	Med-eBase V2 additional user license
383A915	Med-eBase V2 additional product license (Rigel Uni-Therm)



Printers and Scanners

Part No.	Description	Image	Specification
339A923	Bluetooth barcode scanner		Bluetooth scanner for asset detail transfer to Rigel products

Carry Case Solutions

Part No.	Description	Image	Specification
Custom	Rigel Med-eKit Pro		Pelican suitable to hold different Rigel testers (custom to requirements)

At Rigel, we are always developing our range of accessories. If you require an accessory which is not on the list, please contact us.

Rigel Medical, Bracken Hill, South West Industrial Estate,
Peterlee, County Durham, SR8 2SW United Kingdom

Part of  **SEAWARD**
GROUP

Tel: +44 (0) 191 587 8730 **Fax:** +44 (0) 191 586 0227
Email: info@rigelmedical.com **Web:** rigelmedical.com

